

SM-15B Holder Attachment



Stereotaxic micromanipulator SM-15 holds a rod in 5mm O.D. Besides the attachment **SM-15A** can be used with the **SM-15** to hold a rod in 8mm O.D., new attachment **SM-15B** is capable of holding a rod from 3 to 6mm O.D. securely.

The efficient usage of **SM-15B** is to hold the injection holder (4mm O.D.) of IM-300 and other micro-injectors for in-vivo viral and/or dye injection. The existing attachments for **SM-15** require a wrench or a screwdriver to fix a holder. On the other hand, new **SM-15B** does not require any tool. The fixation can be achieved by hand quickly and easily. You might have a doubt on its secure fixation by hand, but it is NOT to worry. The holder within 3-5mm O.D. can be held securely by two steps; clamping and fixing.


The next section shows the procedure of fixing the injection holder in 4mm O.D.

We recommend **SM-15B** for customers who want to fix the holder *quickly, easily, and securely*.

How to fix a holder

The procedure can be achieved by three steps only!

Spring adjustment knob clamps the holder temporarily without a falling, and also becomes useful after its first adjustment due to facilitate the holder replacement.



If you have any questions or need further information, please contact us.