

MOUNTING ADAPTORS FOR STEREOMICROSCOPES

Narishige provides adaptors for stereomicroscopes. These adaptors are either screw-fixed to readily available screw holes on the microscopes or mounted on the illumination column. The screw-fixed type may not be available depending on base stand of the microscope. When you choose an adaptor, please confirm model of the base stand as well as model of the microscope. In this news, we will review the adaptors for stereomicroscopes.

MICROSCOPE	ADAPTOR	AVAILABLE BASE STAND	UNAVAILABLE BASE STAND (SPECIAL-ORDER)
NIKON SMZ645 / 660 / 800/ 1000 / 1500	NN-T NN-T-2	C-PS160 (plain stage) C-DDS (transmitted light stand) C-DSD (transmitted light stand) C-BD (transmitted light stand)	C-DS (transmitted light stand) C-PS / C-PSC (plain stand)
OLYMPUS SZ7 / 61 SZX9 / 10 / 12 / 16	NO-SZX2A	SZX2-ILLT (slim LED transmitted light stand) SZX2-ILLB (transmitted light illuminator) SZX2-ILLD (transmitted light illuminator) SZX2-ILLK (transmitted light illuminator) SZX-ILLD2 (transmitted light illuminator) SZX-ILLB2 (transmitted light illuminator) SZX-ILLB (transmitted light illuminator) SZX-ILLD (transmitted light illuminator) SZX-ILLK (transmitted light illuminator)	SZ2-ILA (transmitted illumination attachment) SZ2-ILST(LED illuminator stand) SZX-ST (standard stand) SZ2-ST (standard stand) SZX-STL (large-size stand)
LEICA MZ series MS5	NL-10-2	MZ series are required illuminator mount type adaptors. The NL-10-2 is available regardless of the base stand. When the MWS-1 series is mounted, the P-6-2 attachment is also required. (Stands requiring the P-6-2 : TLST, TLBF / DF, TLRC, TLRCi) For adaptors for Combi microscopes, please contact us.	
CARL ZEISS Stemi series	NZ-4-2	Stemi series are required illuminator mount type adaptors. The NZ-4-2 is available regardless of the stand. In the case of Discovery and Lumar, the adaptor is special-made. Please contact us for details in this case.	

- * The MMN-1 requires the NR/NR2 with the NO-SZX2A, NL-10-2, NZ-4-2, NN-T-2.
- * Please contact us when you need special-order adaptors.

If you have any questions, please do not hesitate to contact us.